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## Instruction for use MEDISPEED® Liquid medical nitrogen

### Introduction:

- Read all of this leaflet carefully before you start using MEDISPEED.
- Keep this leaflet as you may need to read it again.

### This leaflet gives you information about:

- What is MEDISPEED and Indications/Intended use
- Clinical Contraindications
- Clinical side effects
- How to use Liquid Nitrogen for clinical application
- Safety warnings and advice for the Clinician/Healthcare establishment
- Storage and Transfer of Liquid Nitrogen and Associated Equipment
- Further Information

### What is liquid nitrogen and indications/intended use:

- MEDISPEED is a liquefied gas stored under pressure in purposely designed storage vessels.
- MEDISPEED is intended for use as a cryogen in medical applications due to the cold temperature of the liquid. The boiling point of liquid nitrogen is -196°C.

### Medical Indications/Applications

#### Topical skin conditions (Topical Cryotherapy)

- The most commonplace application for cryotherapies using MEDISPEED is common external skin conditions such as removal of warts, moles, verruca's and/or skin tags. More serious external ailments can also be addressed with this technique, such as solar keratoses and small skin cancers.
- Clinical guidance for Cryotherapy should be used at all times. Clinicians using MEDISPEED for this process should be fully trained and always follow internal procedures.

#### Cancers and other internal conditions (Cryosurgery using Cryoprobes)

- Cryosurgery is rapidly developing as a supplement to traditional cancer therapies. Generally, all tumours that can be reached by cryoprobes are treatable (i.e. solid tumours >1 cm).
- However, the tiny, diffuse metastases that often coincide with cancers are not usually affected by cryotherapy. In addition to cancerous tissues, cryotherapy has also proven a successful therapy to some soft tissue conditions such as plantar fasciitis and fibroma.
- Clinical guidance for Cryosurgery should be used at all times. Clinicians using MEDISPEED for this process should be fully trained in handling MEDISPEED and always follow internal procedures.

#### Cryostorage

- Cryostorage of biological material (cells, tissues, blood, semen, etc.) for subsequent medical procedures and analysis.

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### *Use with other associated devices for purpose of intended use*

- In order to utilise and apply MEDISPEED it is intended for use with other medical devices which are designed to be compatible with MEDISPEED, both in terms of temperature and pressure such as spray applicants (Topical Cryotherapy), probes (Cryosurgery) or freezer storage (Cryostorage).
- Users should always ensure that the aforementioned medical devices they purchase are compatible with MEDISPEED and that they always follow the instructions provided by the manufacturer of these devices.

### **Clinical contraindications:**

#### *Topical Cryotherapy*

- Direct contact application of MEDISPEED such as with a spray style applicator device is not recommended to be used on lesions where the skin is excessively damaged or broken.

### **Clinical side effects:**

- Please be aware that a number of the points detailed in the section may also be a safety hazard for the clinician/user of MEDISPEED.
- Therefore, it is important that MEDISPEED must always be handled by staff trained in and aware of the hazards associated with handling liquid nitrogen.

#### *Topical Cryotherapy*

- Excessive exposure of the application site to the low temperatures of MEDISPEED may cause cold burns and damage/blistering to the surrounding skin. This can be temporary or permanent.
- There may be temporary hypersensitisation and hypopigmentation changes when MEDISPEED is applied directly to skin.

#### *Cryosurgery*

- Side effects from Cryosurgery will depend on the specific tumour being treated. For example, a treatment for skin cancer may leave noticeable scarring, swelling or burning sensations, whereas treatment for prostate cancer may lead to perineal pain and erectile dysfunction, although only in the short term. Tumours located in the cervix may result in cramping, pain or bleeding.

#### *General*

- Transient exposure to very cold gas produces discomfort in breathing and can provoke an asthma attack in susceptible people.
- Low air temperatures arising from the proximity of liquefied atmospheric gases can cause hypothermia and all people at risk should wear warm clothing and minimise time of exposure to the low temperatures.
- However, this is more likely to occur around large quantities of MEDISPEED such as at the storage tanks rather than during clinical use.
- If you experience any of these side effects during use, or if you notice any other side effects not listed in this leaflet, contact Linde Gas Therapeutics Benelux B.V. immediately.

### **How to use Liquid Nitrogen for clinical application:**

- MEDISPEED is intended for use with other medical devices which are designed to be compatible with liquid nitrogen, both in terms of temperature and pressure, such as spray applicants (topical cryotherapy), probes (Cryosurgery) or freezer storage (cryogenic storage).
- Always ensure that equipment used in conjunction with MEDISPEED is compatible with the product.
- Always follow the instructions provided by the manufacturer of these cryogenic applicators/freezer medical devices and storage vessels that are used in conjunction with liquid nitrogen.
- Clinical guidance for Cryotherapy/Cryosurgery should be used at all times. Clinicians using MEDISPEED for this process should be fully trained and follow internal procedures.
- The liquid will vaporise on exposure to atmospheric temperatures and therefore cannot be reused.

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#### Safety warnings and advice for the clinician/healthcare establishment:

- MEDISPEED contains liquid nitrogen which is an extremely cold substance with a boiling point of  $-196^{\circ}\text{C}$  and should be treated with extreme caution at all times. Spills and leaks of MEDISPEED should be avoided.
- Improper handling of MEDISPEED can be fatal and cause personal injury. To prevent this, MEDISPEED must always be handled by staff trained to use liquid nitrogen. Employers should ensure that procedures for the handling of liquid nitrogen are in place and that all staff involved have documented training in these procedures.
- Adequate emergency procedures must be in place in the event of a spill, cold burn or suspected asphyxiation as a result of oxygen deprivation.

#### Embrittlement of materials

- Contact of MEDISPEED with some materials can cause embrittlement which will result in these materials breaking more easily and subsequently pose an increased risk of sharps injuries/cuts when handling the broken material. Where possible these materials should not be used with liquid nitrogen.
- Containers used to store biosamples at cryogenic temperature (either in liquid nitrogen or in the vapours of liquid nitrogen) may become more brittle due the low temperatures and break more easily leading to the possibility of cutting injuries.

#### Other Risks

- Avoid contact of MEDISPEED with the eye(s) by wearing a face mask.
- If contact between MEDISPEED and the eye(s) did occur, rinse the eye(s) immediately with lukewarm water (not exceeding  $37^{\circ}\text{C}$ ) for at least 15 minutes and seek further medical treatment.

#### Storage and transfer of liquid nitrogen and associated equipment:

- As well as the devices used to apply/utilise MEDISPEED in its intended clinical use, storage vessels and dewars for the bulk storage and transport of MEDISPEED to the place of use will be required.
- Again, it is important that these are compatible with liquid nitrogen and comply with all relevant standards and specifications. Linde Gas Therapeutics Benelux B.V. can aid clinicians/healthcare establishments with their choice of storage/transport systems if required.

#### Storage of Liquid Nitrogen

- When choosing a pressurised vessel system for the delivery of and storage of large quantities of MEDISPEED ( $> 200$  litres) make sure the vessels are properly labelled.
- Always follow the instructions of the manufacturer of the pressurised vessel system. Linde Gas Therapeutics Benelux B.V. can aid with the selection, implementation and maintenance of these large storage facilities if required.
- Follow daily checks of the storage system as described within the manufacturer's instructions.
- Always ensure that the level of MEDISPEED within the storage tank is sufficient for your needs. A lack of MEDISPEED within the storage system may result in a delay in treatment options and damage of important samples.
- Ensure that the area where MEDISPEED is stored is equipped with proper ventilation and/or stationary oxygen detection equipment, which will warn in case of the potential danger of oxygen shortage due to the evaporation of liquid nitrogen.
- Specifically, **for the Netherlands**, make sure that liquid nitrogen is stored according to the PGS9 publication "Storage of cryogenic gases of  $0, 125 \text{ m}^3 - 100 \text{ m}^3$ ".

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#### *Maintenance of storage equipment*

- Daily inspection of all MEDISPEED storage equipment must be made by the user to detect any leaks or other malfunctions.
- For more information see the manufacturer's user manual. Service and maintenance of the vessels and accessories should be performed by trained and qualified personnel.
- Contact Linde in case of leakage and unnatural pressure change in the main storage vessel.

#### *Transfer of Liquid Nitrogen between equipment*

- When transferring from the large storage to smaller transport containers/dewars or clinical applicators always follow the instructions of the manufacturer of the equipment.
- Use personal protective equipment and clothing as described in the Gases Care with Cryogenics pamphlet "Safe handling of liquid nitrogen" ([www.linde-healthcare.nl/nl/images/Veiligheid\\_rondom\\_vloeibare\\_stikstof\\_tcm170-596378.pdf](http://www.linde-healthcare.nl/nl/images/Veiligheid_rondom_vloeibare_stikstof_tcm170-596378.pdf)).
- Ensure that people transferring MEDISPEED are fully trained in handling liquid nitrogen and associated equipment.
- Before filling is started, ensure the vessel containing MEDISPEED is stable and unlikely to tip. This is not applicable for the larger storage tanks that will be architecturally fixed in place.

#### *Associated Equipment Use, Storage and maintenance*

- MEDISPEED cannot be stored for any specified period of time in non-pressurized/non insulated transport vessels, dewars or applicators as the liquid will evaporate.
- All containers/applicators used in conjunction with MEDISPEED should be stored in a well-ventilated area to ensure any residual nitrogen escaping to the atmosphere when it evaporates.
- MEDISPEED should not be stored in containers that are not designed to store cryogenic gases. Refer to associated equipment manufacturer's instructions.
- Ice plugs can form in the neck of dewars and can be ejected at high velocity due to pressure build up. Avoid these by ensuring protective caps are used at all times and that dewars are completely emptied before being taken out-of-use or put into storage.
- Safety procedures must be developed and adhered to for the use of dewars including their transportation within and around the premises.
- Special safety procedures are necessary when carrying filled dewars in lifts.
- Only use dewars that are correctly and clearly labelled.
- Always ensure that adequate ventilation is provided in areas where dewars are filled, used or stored.
- Associated equipment must be maintained in accordance with the manufacturer's instructions.
- Please refer to guidelines on "siting of liquid cylinders or vessels in buildings" and "movement of cryogenic vessels in lifts" for further information regarding safe handling of MEDISPEED.

#### *Shelf Life of MEDISPEED*

- Due to heat inleak through the storage vessel walls MEDISPEED will evaporate and will vent to atmosphere on a regular basis. The level of MEDISPEED will be continuously dropping.
- The exact length of time that MEDISPEED will last between deliveries will depend on the insulation and pressure capabilities of the storage vessel and the amount of MEDISPEED that is transferred for use into other containers/dewars/applicators.
- Therefore, it is important that deliveries are scheduled to ensure the amount of MEDISPEED in the storage vessel is maintained for the required usage rate of the liquid.

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Hersteller: Linde Gas Therapeutics Benelux B.V.
Manufacturer: Linde Gas Therapeutics Benelux B.V.